- f. Press **Select**. Small arrows display above and below the setting.
 - g. Press **Up/Down** to select "Printer."
 - h. Press **Select** to save the setting.
 - i. Press **Menu** to close the pop-up and return to the System Menu screen.
 - j. Press **Menu** to return to the monitoring screen.
2. Connect the printer cable to the RS-232 port.
3. When monitoring a patient, an event summary label prints each time the Event Mark button is pressed.

Connecting the Device into a Medical System

Incorporating the device into a medical system requires the integrator to identify, analyze, and evaluate the risks to patient, operators, and third parties. Subsequent changes to the medical system after device integration could introduce new risks and will require additional analysis. Changes to the medical system that must be evaluated include:

- Changing the system configuration
- Adding devices to or disconnecting devices from the system
- Updating or upgrading equipment connected to the system

Issues resulting from user-initiated system changes may include corruption or loss of data.

NOTES:

- Use of a multiple-socket outlet with multiple devices results in a Medical Electrical System.
 - When using the serial port to connect the device to other equipment, follow each device's cleaning instructions.
 - Verify all equipment connected to the device is suitable for the patient's environment.
-

CAUTION: Failure of a network data coupling (serial cable/connectors/wireless connections) will result in loss of data transfer.

Data Output Formats

This device features 5 once-per-second, real-time data output formats (Nonin 1 – Nonin 5). Each data format includes an ASCII header containing model number, time, and date information. In addition, the RS-232 port outputs data through the Dymo printer (Printer).

Formats are selected from the System Menu screen (see "Set Bluetooth and/or RS-232 Data Output Formats" on page 53). Under the Data Output Modes setting, the Bluetooth and RS-232 ports have separate selection options and may use different data output formats.

NOTE: When using SenSmart download software with the monitor, the port used to download data (either Bluetooth or RS-232) must be set to Nonin 1 or Nonin 5 before connecting to the SenSmart software.

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Nonin 1

NOTE: This format is not compatible with all of the X-100M's features.

Baud Rate	57,600
Delimiter	See format in table 13
Line Terminator	CR [0x0D] LF [0x0A]
CRC	CRC-16 CCITT (XMODEM)

Data from the device are sent once per second in the following format:

```
Ch1=XXX Ch2=XXX Ch3=XXX Ch4=XXX 1234&${*|
yyyy-mm-ddThh:mm:ss|rSO2=xxx,xxx,xxx,xxx|HbI=xx.x,xx.x,xx.x,xx.x|
AUC=xxxx,xxx,xxx,xxx|REF=xxx,xxx,xxx,xxx|HI_LIM=xxx,xxx,xxx,xxx|
LOW_LIM=xxx,xxx,xxx,xxx|ALM=xxx,xxx,xxx,xxx|SIG_QUAL_ALM=x,x,x,x|
POD_COMM_ALM=x,x,x,x|SNS_FLT=x,x,x,x|LCD_FLT=x\
LOW_BATT=x\CRIT_BATT=x\BATT_FLT=x\STK_KEY=x\SND_FLT=x\
SND_ERR=x\EXT_MEM_ERR=x\CKSUM=xxxx<CR><LF>
```

NOTE: The 1234&\${* order shall be preserved in all alarm conditions.

Table 14. Nonin 1 Data Output Format

Parameter	Value	Following Delimiter
Ch1=XXX	Channel 1 regional oximeter value. Leading zeros blank; --- if no value available.	space
Ch2=XXX	Channel 2 regional oximeter value. Leading zeros blank; --- if no value available.	space
Ch3=XXX	Channel 3 regional oximeter value. Leading zeros blank; --- if no value available.	space
Ch4=XXX	Channel 4 regional oximeter value. Leading zeros blank; --- if no value available.	space
1234	Patient alarm indication for channel 1,2,3,4. Appears only if a channel's patient alarm is active. If not active, number does not appear (e.g., 14 means channel 1 and 4 patient alarms are active).	none
&	Appears if an equipment alarm is active.	none
\$	Appears if the critical battery state is detected.	none
*	Appears if the event is marked.	
yyyy-mm-ddThh:mm:ss	International date and time format: year, month, day, hour, minutes, seconds.	

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Table 14. Nonin 1 Data Output Format (Continued)

Parameter	Value	Following Delimiter
rSO2=xxx,xxx,xxx,xxx	Regional oximetry values for channels 1,2,3,4 in %. Leading zeros blank; --- if no value available.	
Hbl=xx.x,xx.x,xx.x,xx.x	Hemoglobin index values for channels 1,2,3,4 in grams per deciliter. Leading zeros blank; --.- if no value available.	
AUC=xxxx,xxxx,xxxx,xxxx	Area under curve for channels 1,2,3,4. Leading zeros blank.	
REF=xxx,xxx,xxx,xxx	Reference (low alarm limit) line values for channels 1,2,3,4. Used to calculate AUC. Leading zeros blank.	
HI_LIM=xxx,xxx,xxx,xxx	High limit alarm setting for channels 1,2,3,4. Leading zeros blank. "OFF" if no limit set.	
LOW_LIM=xxx,xxx,xxx,xxx	Low limit alarm setting for channels 1,2,3,4. Leading zeros blank. "OFF" if no limit set.	
ALM=xxx,xxx,xxx,xxx	Active alarm indication for channels 1,2,3,4. Valid values: HI, MAR (marginal), LOW, OFF.	
SIG_QUAL_ALM=x,x,x,x	Signal quality alarm indication for channels 1,2,3,4. 0 = no active alarm. 1 = active alarm.	
POD_COMM_ALM=x,x,x,x	Pod communication alarm indication for channels 1,2,3,4. 0 = no active alarm. 1 = active alarm.	
SNS_FLT=x,x,x,x	Sensor fault indication for channels 1,2,3,4. 0 = no active alarm. 1 = active alarm.	\
LCD_FLT=x	Display fault indicator. 0 = no fault active. 1 = fault active.	\
LOW_BATT=x	Low battery indicator. 0 = no low battery state. 1 = low battery state.	\
CRIT_BATT=x	Critical battery indicator. 0 = no critical battery state. 1 = critical battery state.	\

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Test

- Run algorithm against the string of ASCII characters “123456789”
- Result should be 0x31C3

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Nonin 2

NOTE: This format is not compatible with all of the X-100M's features.

Baud Rate	9,600
Delimiter	Comma [0x2C]
Line Terminator	CR [0x0D] LF [0x0A]
CRC	N/A

Column 1	Column 2	Column 3	Column 4
Current value of Channel 1	Current value of Channel 2	Average of Channel 1 and Channel 2	0

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Nonin 3

NOTE: This format is not compatible with all of the X-100M's features.

Baud Rate	9,600
Delimiter	One or more consecutive spaces [0x20]
Line Terminator	LF [0x0A] CR [0x0D]
CRC	N/A

Version	Date	Time	Channel Name	rSO ₂	Event	Status	Baseline	AUC	UAL	LAL	A	B	C
99.99.99/1/1	mm/dd/yy	hh/mm/ss	These columns repeat per channel										

Continues with:

Sensor ID 1	Sensor ID 2	Sensor ID 3	Sensor ID 4
s/n	s/n	s/n	s/n

Version is 99.99.99/1/1.

Date: month, day, year

Time: hour, minutes, seconds

Channel names:

- L = Channel 1
- R = Channel 2
- S1 = Channel 3
- S2 = Channel 4

rSO₂ readings for missing data outputs as 0.

Event: 0 = no event; 1 = miscellaneous event.

Status values:

- 1 = Pod connected with no sensor (sensor fault)
- 2 = Excessive light indication (used for poor signal quality)
- 4 = Good signal quality, valid rSO₂ readings occurring
- 5 = rSO₂ high alarm
- 6 = rSO₂ low alarm
- 7 = Pod became connected to the system
- 9 = Critical battery alarm
- 11 = Pod is not connected to the system

Baseline and AUC are the current values.

UAL is the upper alarm limit. 0 = OFF.

LAL is the lower alarm limit. 0 = OFF.

A, B, C are 0.

Sensor ID is the sensor image number.

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Event: 0 = no event; 1 = miscellaneous event

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Nonin 5

This data format was designed to be extensible. Future enhancements to the Model X-100M may be included in the data output. As these enhancements become available, new column labels may be added at any position within the data format.

Baud Rate	57,600
Delimiter	Comma [0x2C]
Line Terminator	CR [0x0D] LF [0x0A]
CRC	CRC-16 CCITT (XMODEM)

Date	Time	System Name	Preset Name	Channel Data Columns	Event	Error Code	CRC16
yyyy-mm-dd	hh:mm:ss			The data columns for each connected channel vary depending on the type of sensor connected.		EXX	

Date: year, month, day

Time: hour, minutes, seconds

System name is the name assigned to the system.

Preset name is the name of the preset in use.

Channel data columns: see below for regional oximeter columns, pulse oximeter columns, and unrecognized oximeter columns.

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Test

- Run algorithm against the string of ASCII characters "123456789"
- Result should be 0x31C3

Regional Oximeter Columns

When a regional sensor is attached to the signal processor, the channel data columns are the following:

Channel	Name	rSO ₂	Hbl	AUC	BL	UAL_rSO ₂	LAL_rSO ₂	Status
---------	------	------------------	-----	-----	----	----------------------	----------------------	--------

Channel is the channel number.

Name is the sensor site name.

rSO₂, Hbl, AUC and BL are the current values. Blank if data is missing [dashes on display].

UAL is the current upper alarm limit value.

LAL is the current lower alarm limit value.

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Status is defined by the following and are active high:

Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0
Measurement Accuracy: 0 = Absolute 1 = Trending	Low Alarm: 0 = Auto 1 = Manual	rSO ₂ High	rSO ₂ Low	rSO ₂ Marginal	Sensor Fault	Pod Comm. Lost	Signal Quality

Pulse Oximeter Columns

When a pulse oximetry sensor is attached to the signal processor, the channel data columns are the following

Channel	Name	SpO ₂	PR	UAL_SpO ₂	LAL_SpO ₂	UAL_PR	LAL_PR	Status
---------	------	------------------	----	----------------------	----------------------	--------	--------	--------

Channel is the channel number.

Name is the sensor site name.

SpO₂, PR are the current values. Blank if data is missing [dashes on display].

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Status is defined by the following and are active high:

Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0
Low Perfusion	PR High	PR Low	SpO ₂ High	SpO ₂ Low	Sensor Fault	Pod Comm. Lost	Signal Quality

Unrecognized Oximeter Columns

When there is no sensor attached to the signal processor, the channel data columns are the following

Channel	Name	Status
---------	------	--------

Channel is the channel number.

Name is the sensor site name.

Status is defined by the following and are active high:

Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0
Unused	Unused	Unused	Unused	Unused	Unrecognized Sensor	Pod Comm. Lost	Unused

Printer

This data output format was designed to work with the Dymo LabelWriter SE450 Label Printer. See "RS-232 Connection to Printer" on page 63 for more information.

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